


TOTAL KNEE SYSTEM

#1600-13, -14, -16
SET-UP INSTRUCTIONS

Thank you for purchasing the Sawbones TKS - Total Knee System. This versatile knee holder was designed for demonstrating and training total knee arthroplasty on medium and large knee models. The TKS allows tibial alignment systems to register on the foot, which features anatomical lateral and medial malleoli for digitization of surgical navigation. Purchasing options and setup instructions are outlined below.



- ①
travel bag
- ②
foot clamp
- ③
pelvis clamp
- ④
bone clamp
- ⑤
acetabular holder
(medium or large)
- ⑥
soft tissue
foot
- ⑦
full length knees
(sold separately)

PURCHASE OPTIONS

TKS
#1600-13 ① ② ③ ④ ⑤ ⑥
#1600-14 ① ② ③ ④ ⑤ ⑥
#1600-16 ① ② ③ ④ ⑤ ⑥

COMES WITH
LARGE ONLY
MEDIUM ONLY
BOTH MEDIUM & LARGE

RELATED PARTS ②

Additional foot clamps can be purchased individually as #1600-17.

FULL LENGTH KNEES (SOLD SEPARATELY) ⑦

#1145 & #1145-1 (large)
use with #1600-13, #1600-16
#1179 & #1179-1 (medium)
use with #1600-14, #1600-16

SET-UP INSTRUCTIONS

