

MODEL SPECIFICATIONS

OVERVIEW

MODEL NAME	Ultrasound Neonatal Internal Jugular and Carotid Cannulation Trainer
SKUS	500-303-30 (clear tone) 500-303-32 (medium tone)
MANUFACTURER	Sawbones
BRAND	RediStik



CANNULA/DILATOR

JUGULAR VEIN	≤ 15 French
CAROTID ARTERY	≤ 9 French
RESONANCE™ TISSUE	≤ 23 French



DURABILITY

JUGULAR VEIN	5-10 Cannulations
CAROTID ARTERY	5-10 Cannulations
RESONANCE™ TISSUE	~ 50 Cannulations

KEY FEATURES

- Supports the Seldinger technique (wire-guided, ultrasound-guided) for neonatal HD, CVC, and ECMO catheter placement
- Self-healing Resonance™ tissue shows little to no leaking with proper care
- Vessels are independently replaceable using forceps and lubricant (Ultrasound Neonatal Internal Jugular Vessel Replacement, 599-50-1; Ultrasound Neonatal Carotid Vessel Replacement, 599-50-2); tissue replaceable with the Ultrasound Neonatal Internal Jugular and Carotid Cannulation Replacement, Clear Tone (525-303-30), Medium Tone (525-303-32)
- Can be used standalone or under ultrasound guidance



MADE WITH
resonance™
 A SAWBONES® TISSUE