

MODEL SPECIFICATIONS

OVERVIEW

MODEL NAME	Ultrasound Branching 3-Vein Access Trainer with Artery
SKUS	500-502-30 (clear tone) 500-502-32 (medium tone)
MANUFACTURER	Sawbones
BRAND	RediStik
NEEDLE/CATHETER	≤ 22 AWG
DURABILITY	over 1,000 sticks



KEY FEATURES

- ▶ Ideal for introductory training and instructor-led demonstration
- ▶ Features a 5 mm bifurcated vessel with a 3 mm branch
- ▶ Includes a 3 mm vein with valve and non-collapsing 4 mm artery
- ▶ Supports flashback confirmation and simulated blood draw
- ▶ Can be used with or without ultrasound
- ▶ Allows for flushing and infusion of fluids
- ▶ Compatible with alcohol, CHG, and dressings for simulated skin prep and device securement*
- ▶ Kit includes setup and transport accessories
- ▶ Resonance™ material is heat resistant and can withstand temperatures up to 150° F without risk of deformation

*all devices must adhere to warranty policy and not be left in the model for extended periods of time



MADE WITH
resonance™
 A SAWBONES® TISSUE