



MADE WITH
resonance™
A SAWBONES® TISSUE

MODEL SPECIFICATIONS

OVERVIEW

MODEL NAME	Clear Heart Visualization Module
SKUS	501-301
MANUFACTURER	Sawbones
BRAND	RediStik

KEY FEATURES

- ▶ Rigid, fully transparent construction for unobstructed visualization
- ▶ Features all four heart chambers
- ▶ Locks onto the Ultrasound Internal Jugular and Carotid Cannulation Module (500-301-30/32) for single- or dual-site ECMO simulation
- ▶ Also compatible with the 500-702 series femoral access model (COMING SOON)
- ▶ Manifolds allow for easy priming in any configuration
- ▶ Dedicated air-vent port keeps the chambers clear for visual inspection
- ▶ Ideal for instructor-led demonstration and visual confirmation of catheter placement
- ▶ Supports full-circuit ECMO cannulation